

Lipid Peroxides (LPO) Kit for biological samples




- For the quantitative determination of high levels of lipid peroxides in serum, plasma, and other biological samples
- Specifically quantitates lipid peroxides. Does not measure metabolites of lipid peroxides
- Applicable to clinical chemistry analyzers
- Manual or automated applications
- 1 point calibration, 15 minute assay
- Number of tests depends on volume of cuvette:

3 mL cuvette = 55 tests

1 mL cuvette = 165 tests

Microtiter well = 550 tests



Assay	Components	Assay Range	Wavelength	Catalog No.
LPO-CC	Reagent 1: Enzyme, 55 mL Reagent 2: Chromogen, 110 mL Calibrator: 50 nmol/mL, 10 mL	2 - 300 nmol/mL	675 nm	CC-004

For *In Vitro* Research Use Only. Not for use in diagnostic procedures in the